



Electrospinning Accessories



Electrospinning accessories will improve your experiment through a safer, more versatile and mobile electrospinning. The electrospinning chamber provides an enclosed area to protect your electrospinning from environmental variables such as air flow, and also improve electrospinning safety by providing a protective layer between you and the high voltage. The equipment trolley enable a mobile electrospinning setup with the versatility to place the electrospinning components according to your own arrangement. Additional accessories such as the coaxial nozzle enables you to produce core-shell structure of nanofibers, while the blunt needle is a hard-to-miss accessory of many standard electrospinning setups.




Chamber

- ◇ Electrospinning chamber for basic series.
- ◇ Coloured chamber equipped with LED lighting for easy observation.
- ◇ Exhaust fan and duct connection ready.
- ◇ Removable transparent front lid and easy access to electrospinning components

Trolley

- ◇ Trolley for basic series.
- ◇ Mobile trolley with castor wheel.
- ◇ Neat arrangement of components.
- ◇ Dimension(mm): 900(L) x 700(W) x 850(H)



Attachable vertical plat collector

- ◇ Stainless steel plate collector designed to attached to drum collector.
- ◇ Easy interchange with drum collector during electrospinning.
- ◇ Size: 31cm x 25cm

Stand alone vertical plate collector

- ◇ Stainless steel plate collector completed with base and grounding wire.
- ◇ Size: 25cm x 20cm

** Product specifications are subject to change without prior notice due to engineering innovations*